


DC and SCP Cycle Training

DC and SCP Cycle Training

This design is based upon the “Tubingen” approach of training SCP up and down alternately.





You must have the Session Wizard to do this. MINI-Q owners have this. Others have to get a license for the Session Wizard.

DC and SCP Cycle Training

The feedback changes every 8 seconds. For 8 seconds, you hear a high sound when the SCP is rising. Then for 8 seconds you hear a low sound when the SCP is falling. This repeats 120 times.


Flash games and animations and DVD's, CD's, etc are enabled. The rewards come when either task is being successfully done, so as the trainee alternates to increase and decrease SCP, they build up rewards in the game.



Puzzle and Space Race will respond to the points counter 1.

DC and SCP Cycle Training


When starting the session (Run the Next Session / GO) First you will see this:



Then this:




Just "click through" these. Then press OK when you see good signals on the screen below.



DC and SCP Cycle Training

This design uses “Full Atlantis” mode combined with the Event Wizard and the Session Wizard.


Flash Games & multimedia are controlled by Event #5



Points are on events 3 (up) and 4 (down) – can add refractory period or sustained reward if desired
MIDI is on events 14 (up) and 15 (down) – can adjust sounds voices, pitch, modulation, etc.

DC and SCP Cycle Training

The events are listed below:





The screenshot shows a window titled "BMrEventGUI30" with a yellow warning icon in the top-left corner. The window contains a list of eight events, each with a summary of its logic and its current status. The events are listed as follows:

Event Summary	Current Status
Summary for Event 1: IF: EQN: x=E11A-E11B; // SCP IS GREATER THAN Value: 0.0 SUSTAIN: 100ms THEN: Do Nothing	EVENT 1 IS CURRENTLY: ENABLED
Summary for Event 2: IF: EQN: x=DC1; IS GREATER THAN Value: 0.0 SUSTAIN: 100ms THEN: Do Nothing	EVENT 2 IS CURRENTLY: ENABLED
Summary for Event 3: IF: EQN: x=E14F; IS GREATER THAN Value: 0.0 REFRACTORY: 500ms THEN: Award Point (Counter 1)	EVENT 3 IS CURRENTLY: ENABLED
Summary for Event 4: IF: EQN: x=E15F; IS GREATER THAN Value: 0.0 REFRACTORY: 500ms THEN: Award Point (Counter 2)	EVENT 4 IS CURRENTLY: ENABLED
Summary for Event 5: IF: EQN: x=E16F; IS GREATER THAN Value: 0.0 SUSTAIN: 100ms THEN: Do Nothing	EVENT 5 IS CURRENTLY: ENABLED
Summary for Event 6: IF: Value: 0.0 IS GREATER THAN Value: 0.0 SUSTAIN: 100ms THEN: Do Nothing	EVENT 6 IS CURRENTLY: DISABLED
Summary for Event 7: IF: Value: 0.0 IS GREATER THAN Value: 0.0 SUSTAIN: 100ms THEN: Do Nothing	EVENT 7 IS CURRENTLY: DISABLED
Summary for Event 8: IF: Value: 0.0 IS GREATER THAN Value: 0.0 SUSTAIN: 100ms THEN: Do Nothing	EVENT 8 IS CURRENTLY: DISABLED

An "OK" button is located at the bottom center of the window.

DC and SCP Cycle Training

BMrEventGUI30 

 Summary for Event 9:
IF: EQN: x=E11A-E11B; IS GREATER THAN Value: 0.0 SUSTAIN: 100ms
THEN: Do Nothing
EVENT 9 IS CURRENTLY: ENABLED

Summary for Event 10:
IF: EQN: x=E11A-E11B; IS LESS THAN Value: 0.0 SUSTAIN: 100ms
THEN: Do Nothing
EVENT 10 IS CURRENTLY: ENABLED

Summary for Event 11:
IF: EQN: x=C1DCP; IS LESS THAN EQN: x=C1DCP;SUSTAIN: 100ms
THEN: Do Nothing
EVENT 11 IS CURRENTLY: ENABLED

Summary for Event 12:
IF: EQN: x=FEEDBACK; IS GREATER THAN Value: 0.0 SUSTAIN: 100ms
THEN: Do Nothing
EVENT 12 IS CURRENTLY: ENABLED

Summary for Event 13:
IF: EQN: x=FEEDBACK; IS LESS THAN Value: 0.0 SUSTAIN: 100ms
THEN: Do Nothing
EVENT 13 IS CURRENTLY: ENABLED


Summary for Event 14:
IF: EQN: x=E9F+E12F; IS GREATER THAN Value: 1.0
THEN: Play MIDI Sound MODE: 2 NOTE: 44 INSTR: 48 Strings
STYLE: Sustained MODULATION: Amplitude LOUDNESS: Level: 100
LOUDNESS CHANGE RATE: 3 PITCH CHANGE RATE: 3
KEY: A MODE: Chromatic CHORD: 1 Note
EVENT 14 IS CURRENTLY: ENABLED

Summary for Event 15:
IF: EQN: x=E10F+E13F; IS GREATER THAN Value: 1.0
THEN: Play MIDI Sound MODE: 2 NOTE: 37 INSTR: 16 Organ 1
STYLE: Sustained MODULATION: Amplitude LOUDNESS: Level: 100
LOUDNESS CHANGE RATE: 3 PITCH CHANGE RATE: 3
KEY: A MODE: Chromatic CHORD: 1 Note
EVENT 15 IS CURRENTLY: ENABLED

Summary for Event 16:
IF: EQN: x=E14F+E15F; IS GREATER THAN Value: 0.0 SUSTAIN: 100ms
THEN: Do Nothing
EVENT 16 IS CURRENTLY: ENABLED


DC and SCP Cycle Training

SCP data can be reviewed using BReview: (Review Session Results)
Turn off all the “Channels” and select Event 1



DC and SCP Cycle Training

This is not real data, but it shows how the display works. We still have to work on scaling the graph for negative numbers.



Note that during this training, the BrainMaster continues to acquire all EEG data including alpha, beta, gamma amplitudes, and you can even use z-scores during this as well.

This program can be easily increased to 2 or 4 channels by modifying some of the Events. The basic design would not change.

DC and SCP Cycle Training

The Session Wizard program looks like the following:

BrainMaster Session Wizard Control File

VERSION: 2

NAME: DC and SCP Feedback

COM: DC and SCP Feedback

NCHANS: 1

BEGCYCLE: 1

FEEDBACK: 1

DUTY: 100

PERIOD: 8

FEEDBACK: -1

DUTY: 100

PERIOD: 8

ENDCYCLE: 1

DO:

DOCYCLE: 1

DOCYCLE: 1

DOCYCLE: 1

(total of 120 times)