



# BRAINMASTER MULTI MEDIA PLAYER WITH BUILT-IN DVD PLAYER

**THE BRAINMASTER MULTIMEDIA PLAYER SOFTWARE** is an optional add-on tool box that also provides advanced animations and sounds for neurofeedback training. It is integrated into the BrainMaster 2.5SE and 3.0, Discovery 1.4, and BrainAvatar™ 4.0 software and is the standard option for multimedia and animations. It can play animation files (avi, mpg, mpeg, or wmv format), and will automatically adjust the animation to suit whatever protocol is being used. The system automatically detects your protocol settings, and makes the selected animation run faster or slower, depending on the training results.

Up to 50 frames per second can be shown, with moment-by-moment control of the feedback display. In addition, any sounds built into the animation will also be played, and will also stop and start along with the animation, for a full audio and visual effect. Dozens of animations are included with the basic MMP software license, and you can add your own at any time. There are many sources of animations, both licensed and as freeware, to choose from. Any files can be simply copied into the MMP directory on your PC, and will be immediately usable by the MMP software.

In addition, the MMP player can play sound files and CDs, using a unique “amplitude modulated” method. Sound files (wav, mp3, mp4, or midi) can be selected and played, regardless of the animation files being used. The system can also be used for sound-only feedback, using saved sound files. The MMP player can also play CDs directly. Just insert a CD; select the track(s), and the MMP player will provide EEG-controlled sound feedback.



All DVDs must be obtained from an authorized distributor. Rio ©2011 Twentieth Century Fox/Blue Sky Studios. Carlos Saldanha, Director

The MMP player provides separate controls for the built-in sounds in your animation files, and for the separate files or CDs that you use. This operates in addition to the standard sounds produced by the BrainMaster software. Thus, you can have very complex and engaging sounds, creating layers of feedback for your trainees. If you want to view animations that are on DVDs, you can do so immediately with the built in DVD player.

**COMPATIBILITY:** Windows XP, Vista, or 7 required. BrainMaster MMP requires you to have Microsoft Direct X 9.0c or above installed on your PC (A DX9.0c installation file is included with the MMP software). MMP is compatible with Gforce and similar sound-driven graphics extension programs, as explained in the MMP manual. Alpha/Theta training software and sound program also included!

BrainMaster Multimedia/DVD Player (533-401)	\$450.00 (\$250.00)
Also available for Discovery NeuroGuide MMP (533-416)	\$450
Complete Multimedia Package for Brainmaster and NeuroGuide (533-610)	\$900 (\$650)

A COMPATIBLE THIRD PARTY DVD DECODER IS REQUIRED for proper audio and video playback with the BMR DVD player. Please see <http://www.brainm.com/kb/entry/119/> for specific requirements and instructions for Windows XP, and [www.brainm.com/kb/entry/121/](http://www.brainm.com/kb/entry/121/) for Windows Vista and 7.



**195 Willis Street  
Bedford OH 44146  
440.232.6000 ext 200 or 201  
www.brainmaster.com  
sales@brainmaster.com**

