

BRAINDRYVR CASCADE



Use on any EEG Feedback device.

“BRAINDRYVR” CASCADE (395-002): Offer clients an experience they won’t soon forget. Includes power adapter and signal cable. The Braindryvr Cascade is a stand-alone unit designed to be used with any EEG feedback system. The unit is a sound operated gate that permits the control of electrically operated devices contingent on any EEG event that the feedback system allows. Because the units are “cascadable” (multiple units can be controlled by the same EEG event) the therapist can design highly sophisticated multiple contingency feedback protocols. The supersensitive contingency control permits qualitative differences in EEG events to have different contingencies. Other features offer the therapist many options for unique protocol design. The sound trigger coming from the EEG feedback system has a direct output, which can be used to provide the client with direct EEG feedback when a threshold is satisfied but some other stimulus condition if the threshold is not satisfied, or as a trigger for additional cascaded brain driving units. The delay function of the Braindryvr Cascade permits control of stimulus duration so that the time the stimulus remains on after the triggering event can be controlled by the therapist. Perfect reward for kids!

W/ adapter \$475.00 (\$395.00) W/o adapter deduct \$15.00.

BrainMaster Technologies, Inc.
24490 Broadway Avenue
Oakwood Village, Ohio 44146
Corporate Offices (440) 232-6000 Fax (440) 232-7171
info@brainm.com